



# Carrot C-COM Pro

## WIRELESS HEADSET SYSTEM



**Base System Includes:** (4) HD Wireless Headsets, (1) Wired Headset, (1) Battery Charger, (8) Batteries, and a **2 Year Warranty**

**Benefit from greater clinical workflow efficiency with flawless audio fidelity and precise communication among members of the clinical team.**

*“The communication system improves the entire collaboration of the procedure. I feel that by having everyone on the team listening on the headsets, people are more engaged and even anticipate my needs. The collaboration of the case is vastly improved.”*

## Features

- \*  
Whisper-sensitive individual microphones
- \*  
Flawless audio fidelity
- \*  
Battery life up to 20 hours
- \*  
300 foot range
- \*  
Expandable up to 15 headsets

## C-Com in Action

- London Health Science Centre
- Trillium Health Centre
- Hospital of the University of Pennsylvania
- Albert Einstein Medical Center
- UMass Memorial Medical Center
- Yale New Haven Hospital
- Brigham and Women’s Hospital
- Penn Presbyterian Medical Center

**Plus 100s more ...**

# MedTach

1-289-644-4985  
 info@medtach.com  
 www.medtach.com

# C-Com Pro Specifications



## PERFECT AUDIO. PRECISE

### ➤ System:

Frequency .....	2.4-2.4835 GHz FHSS
Concurrent Talk Paths .....	5
Total Listen paths .....	15
Speaker Output .....	RCA jack
Dimensions .....	2U, W20.4" x H5.0" x D15.9"
Power Requirement .....	AC 100-240V 50/60 Hz



### ➤ Headsets:

Battery Life .....	Up to 15 hours
Weight .....	3.7 oz.
Microphone .....	Directional Noise Cancelling
Battery .....	Removable rechargeable Lithium Ion (2 per headset)
Low Battery Indication .....	LED Blink & Voice Prompt
Approximate Operating Range .....	Up to 300 feet
Ear piece Covers .....	200 Disposable

### ➤ Battery Charging Station:

Charging Bays .....	4 active, 4 reserve
Power .....	Max. 500W (500VA)
Dimensions .....	W5.5" x H3.0" x D7.0"
Weight .....	2.0 lbs.
Operating Temperature .....	5-35° C

***MedTach***

Bringing Innovations to Canadian Clinicians and Researchers

Phone: 1-289-644-4985

Email: [info@medtach.com](mailto:info@medtach.com)

[www.medtach.com](http://www.medtach.com)