


# Cap Sizing Chart

## **MedTach Inc.**

4145 North Service Road, Suite 200,  
Burlington, Ontario L7L 6A3  
Phone 289-644-4985  
[info@medtach.com](mailto:info@medtach.com)



For an accurate head measure, place a tape measure across the forehead and measure around the full circumference of the head. Keep the tape snug for accurate results.



Hat Size	Size	Inches	Cm
6 3/8	X-Small	20.7	53
6 7/8	Small	21.0	54
6 7/8	Small	21.4	55
7	Medium	21.8	56
7 1/8	Medium	22.2	57
7 1/4	Large	22.6	58
7 3/8	Large	23	59
7 1/2	X-Large	23.4	60
7 5/8	X-Large	23.8	61
7 3/4	XX-Large	24.2	62
7 7/8	XX-Large	24.6	63
8	XXX-Large	25.0	64