

## Who we are

MedTach is a Canadian medical device distribution company dedicated to delivering the highest level of customer service to physicians, nursing, technicians, researchers, and ultimately the patients whose treatment regime includes products supplied by the professionals at MedTach.

We focus on delivering innovative technologies including the advanced auscultation solutions provided by the Eko Core Stethoscope.



Cardiology



Primary Care



Education



Telemedicine

## Contact Us

If you have questions or comments about the Eko Core Stethoscope please contact MedTach at:

Phone: 1-289-644-4985

Email: [info@medtach.com](mailto:info@medtach.com)


Web: [www.medtach.com](http://www.medtach.com)



4145 North Service Rd.,  
Suite 200,  
Burlington, Ontario, Canada  
L7L 6A3

# Eko Core Smart Stethoscope

The most versatile digital stethoscope



Bringing Innovations to Canadian Clinicians and Researchers

[www.medtach.com](http://www.medtach.com)

# Eko Core Smart Stethoscope

- ✓ Switchable between 'regular' analog and amplified modes
- ✓ 40X amplification in digital mode
- ✓ Crystal-clear digital audio with white noise reduction
- ✓ Works with heart, lung & other body sounds with visualizations
- ✓ Rechargeable battery with 21-days use average between charges
- ✓ Bluetooth smartphone/tablet connectivity
- ✓ Compatible with your current stethoscope
- ✓ Use as telemedicine stethoscope and stream audio in real-time



Visualize, record, and share sounds wirelessly from your Eko Core to your **Apple** or **Android** mobile device in real time with our FREE app.

40x

Amplification



Noise Filter



Wireless Connection

## The most versatile stethoscope

### Toggle between Analog & Digital




Eko Core allows you to quickly switch between analog (off) and amplified (on) modes. Keep the analog sound but gain all the functionality of digital. Eko is the only stethoscope that allows the ability to switch between analog and digital modes.

### Intuitive Software



We've built the world's first Bluetooth-connected auscultation app for mobile devices. The Eko app streams audio wirelessly from your Eko digital stethoscope to help you confidently and quickly assess patient's heart, lung and body sounds. Visualize, record, share, and save critical sounds from your stethoscope.

### Use Eko Core with your current stethoscope



Select either the Eko Core complete electronic cardiology-grade stethoscope, or the Eko Core Attachment that connects to (and digitizes) the analog stethoscope you already own.